




**Do you have a Validated Biomarker for this
Compound?**

Kirk Newland
Raymond Farnen
Rafiq Islam

Do you have a validated biomarker for this compound? Questions that need to be answered

- Is this for a small molecule or a large molecule?
- What is the species? Is this the right biomarker for this species?
- What is the intra-subject and inter-subject variability?
- **Is this data going to be used for a PK/PD plot or for a statistical comparison?**
- Is this an endogenous biomarker? Will this biomarker be present in the matrix? Do different disease states affect the concentration of this biomarker?
- **Is this biomarker present in the environment?**
- Is reference material available for this biomarker?
- **How much are you are willing to spend for this assay?**

Presentation Outline

- Role of Biomarkers in MRTTP development
 - Fit-for-purpose Bioanalytical Method Validation
 - The continuum of biomarker validation
 - Examples where improvements to the method can result in decreasing the clinical study size
 - Sample Collection
 - Assay sensitivity
 - Assay precision
- 
- Conclusion: Do you have a validated biomarker for this assay that will reduce clinical costs?
 - The ability of the biomarker assay to affect the clinical study is exciting!

Modified Risk Tobacco Product (MRTP) - FDA Guidance

- Exposure Modification Order (Section 911(g)(2)) requires
 - Scientific evidence of substantial overall reductions in exposure to the harmful substance(s)
 - The available scientific evidence demonstrates that a measurable and substantial reduction in morbidity/mortality among individual users is reasonably likely in subsequent studies
 - The order would be appropriate to promote the public health


Challenges to Establishing Reduced Harm using section 911(g)(2)

- MRTP development requires faster and more controllable methods to assess smoking related diseases other than epidemiological studies
- Biomarkers offer an alternative and cost-effective approach for evaluation of potential harm reduction from MRTPs during product development
- The profiles of biomarkers may be used to understand biological events from smoke inhalation to disease manifestation
- Endpoint is the biological effect in response to smoking as opposed to disease manifestation

How SMART is your Biomarker?

- As suggested by Shehabi and colleagues a SMART biomarker is:
 - Sensitive (and selective)
 - Measurable (with high precision)
 - Available (affordable and safely attainable)
 - Responsive (and reproducible)
 - Timely (to allow decision making)
- And to be SMARTER one adds Evaluate (validate) and Re-evaluate (revalidate)


Fit for Purpose Biomarker Validation





When biomarker data will be used to support a regulatory action, such as the pivotal determination of safety and/or effectiveness or to support labeled dosing instructions, the assay should be fully validated in a bioanalytical lab.
- FDA BMV Guidance 2013

Fit for Purpose Biomarker Validation

What is the purpose for each assay?



Is the value normal or abnormal?



Must use the right assay for the right purpose:

- Example (CDER)
 - Company was developing a diabetes drug
 - Conducted a clinical study – GOAL: demonstrate a statistical decrease in blood glucose levels
 - Used a diagnostic clinical chemistry assay because it was cheaper
 - FDA rejected the study
 - Study was rerun and samples were analyzed using a bioanalytical biomarker assay
 - Study was accepted by the FDA
- Lesson learned for section 911(g)(2) filing?

iPF_{2a}-III Oxidative Stress Biomarker

- Identify the “validated method” you would like to use to evaluate tobacco smoker exposure risk as a primary clinical end point.

Study	Detection Method	Concentration of Healthy Normals (ng/mg creatinine)
Sircar et al.	LC-MS	0.36 (0.13)
Ohashi. et al.	LC-MS	0.43 (0.06)
Mori et al.	GC-MS	1.14 (0.02)
Lee et al.	GC-MS	0.50 (0.40)
Tsikas et al.	GC-MS	0.29 (0.04)
Liang et al.	LC-MS/MS	0.25 (0.15)

Bioanalytical Tobacco Assay Method Validation

- Bioanalytical Guidance – assay must be:
 - Appropriately collected ✓
 - Selective
 - Sensitive ✓
 - Precise ✓
 - Stable
 - Sample collection and handling
 - Freezer (-20°C or -80°C)
 - Freeze/Thaw
 - Benchtop
 - Pre-extraction
 - Post-extraction
 - DOCUMENTED !

Method Validation – Sample Collection

- Stability to light
 - Mercapturic Acids and some Tobacco Specific Nitrosamines are light sensitive in urine.
 - Must collect and transport the urine samples in light impenetrable tubes
- Minimize environmental contamination:
 - o-toluidine, an aromatic amine, was shown to leach from certain urine collection containers (4X > LLOQ) and other types of plastic ware.
 - Nicotine is ubiquitous (tube contamination 20X > LLOQ) - Methanol rinsed tubes, pipettes, etc. eliminates contamination.



Method Validation – Sample Collection


COST SAVINGS

- If you don't collect the samples properly to stabilize them and to eliminate contamination then don't bother to run the study.
- Proper collection and sample handling requires the preparation of a detailed sample handling manual.

Method Validation – Sensitivity

EXAMPLES: NNN & 3-OH-B[a]P


- Old LLOQ (most sensitive)
 - NNN = 0.70 pg/ml
 - 3-OH-B[a]P = 50 fg/ml
- Old assay results (% of clinical samples that were < LLOQ):
 - NNN: >20%
 - 3-OH-B[a]P: > 40%
- Solution: purchased an AB Sciex 6500 LC-MS/MS and developed a new assay:
 - NNN LLOQ = 0.20 pg/ml
 - 3-OH B[a]P = 25 fg/ml
 - Dramatic decrease # samples < LLOQ



If a sample is < LLOQ then that subject cannot be included in the statistics


Method Validation – Sensitivity

COST SAVINGS: Able to dose fewer participants




REVIEW: Fit for Purpose Biomarker Validation



What is the purpose for each assay?



Is the value
normal or
abnormal?




Precision - Assay variability comparison



The design of bioanalytical assays using multiple standards produces constant variability over a large concentration range.


Method Validation – Precision

COST SAVINGS: Fewer Participants



Method Validation: Sensitivity + Precision

COST SAVINGS: Almost additive



Biomarker Evolution

Time



Biomarker Discovery

Validated Biomarker

Biomarker optimized for clinical program

- Testing 6 -12 biomarkers / assay
- Clinical Chemistry or Luminex for multiplexing
- High variability
- Sensitivity has not been determined
- Selectivity is not validated

- Testing 2 -5 biomarkers / assay
- OK variability
- OK sensitivity
- Selectivity is confirmed
- Meets regulatory guidelines

- Testing 1- 4 biomarkers / assay
- Low variability
- Sensitivity (linear range) has been optimized for MRTTP
- Customized for program / product

Conclusion

- Do you have validated biomarker for this assay?
 - Multiple answers to that question
 - Not all validated biomarker assays are equal
- The performance of biomarker assays for statistical comparisons for MRTP development can have a direct impact upon the number of subjects required for the clinical program.

References

- “Biomarkers in Acute Lung Injury – Marking Forward Progress” Nicolas Barnett, et al. Crit Care Clin **2011**
- “Isoprostane Measurement in Plasma and Urine by Liquid Chromatography – Mass Spectrometry with One-Step Sample Preparation” Debajit Sircar, et al. Clinical Chemistry **2007**
- “F2-Isoprostanes as Novel Biomarkers for Type 2 Diabetes: a Review” Subremanian Kaviarasan et al. J Clin Biochem Nutr. **2009**