



Global List of Validated Bioanalytical Assays

For further information on our validated assays, please contact your
business development representative or via info@celerion.com.

Assays Under Development

Compound Name	Technique	Species/Matrix
1-AMN; 2-AMN; 3-ABP; 4-ABP	LC-MS/MS	Human Urine
Aprepitant	LC-MS/MS	Human Plasma
Fluoxetine; Duloxetine; Paroxetine; Sertraline	LC-MS/MS	Human Plasma
Nornicotine; Norcotinine; Nicotine-N'-oxide; Cotinine-N-oxide; 4-Hydroxy-4-(3-pyridyl)-butanoic acid	LC-MS/MS	Human Urine
Octreotide	LC-MS/MS	Human Plasma
o-Toluidine	LC-MS/MS	Human Urine
Propylene Glycol	LC-MS/MS	Human Plasma

Validated Assay List

Compound Name	Technique	Species/Matrix	Range	Anti-coag.
1-Hydroxypyrene (1-OHP)	LC-MS/MS	Human Urine	10-2000 pg/mL	N/A
1; 2-Dihydroxybutyl Mercapturic Acid (DHBMA)	LC-MS/MS	Human Urine	10-1000 ng/mL	N/A
2-Pralidoxime (2-PAM)	LC-MS/MS	Human Plasma	0.1-20 ng/mL	Heparin
3-Hydroxypropylmercapturic Acid (3-HPMA)	LC-MS/MS	Human Urine	35-5000 ng/mL	N/A
5-Fluorocytosine	LC-MS/MS	Human Plasma	0.5-100 ug/mL	EDTA
5-Fluorouracil	LC-MS/MS	Human Plasma	25-5000 ng/mL	EDTA
6-Beta Naltrexol	LC-MS/MS	Human Plasma	0.5-25 pg/mL	EDTA
6-Mercaptopurine (6-MP)	LC-MS/MS	Human Plasma	0.25-25 ng/mL	Heparin
8-iso-Prostaglandin F _{2α} (Type III)	LC-MS/MS	Human Urine	25-1500 pg/mL	N/A
8-iso-Prostaglandin F _{2α} (Type IV)	LC-MS/MS	Human Urine	25-6000 pg/mL	N/A
11-Dehydro Thromboxane B ₂	LC-MS/MS	Human Urine	5-500 pg/mL	N/A
17-Deacetyl norgestimate; Levonorgestrel	LC-MS/MS	Human Plasma	20-2500 pg/mL; 80-3500 pg/mL	EDTA
Abacavir; Zidovudine	LC-MS/MS	Human Plasma	25-4000 ng/mL; 25-4000 ng/mL	Heparin
Abacavir; Zidovudine; Lamivudine	LC-MS/MS	Human Plasma	30-9000 ng/mL; 10-3000 ng/mL; 10-3000 ng/mL	EDTA
Acetamidophenol	LC-MS/MS	Human Plasma	0.5-25 ug/mL	Heparin, EDTA
Acetaminophen	LC-MS/MS	Human Plasma	50-10000 ng/mL	EDTA, Lithium Heparin
Acetyl Salicylic Acid (ASA)	LC-MS/MS	Human Plasma	0.1-20 ug/mL	EDTA
Albuterol	LC-MS/MS	Human Plasma	10-2000 pg/mL	K ₂ EDTA
Alendronate	LC-MS/MS	Human Urine	1-200 ng/mL	N/A
Alendronate	LC-MS/MS	Human Serum	0.05-20 ng/mL	N/A
Alendronate	LC-MS/MS	Human Plasma	50-20000 pg/mL	EDTA
Allopurinol; Oxypurinol	LC-MS/MS	Human Plasma	0.1-20 ug/mL	EDTA
Alprazolam; Alpha-Hydroxyalprazolam	LC-MS/MS	Human Plasma	0.10-15 ng/mL; 0.05-7.5 ng/mL	K ₃ EDTA
Amiodarone; Desethylamiodarone	LC-MS/MS	Human Plasma	10-10000 ng/mL; 0.5-200 ng/mL	EDTA
Amlodipine	LC-MS/MS	Human Serum	50-10000 pg/mL	N/A
Amlodipine	LC-MS/MS	Human Plasma	50-10000 pg/mL	EDTA

Validated Assay List

Compound Name	Technique	Species/Matrix	Range	Anti-coag.
Amoxicillin	LC-MS/MS	Human Plasma	0.1-30 ug/mL	EDTA
Amphetamine D/L	LC-MS/MS	Human Plasma	0.5-50 ng/mL	EDTA
Anastrozole	LC-MS/MS	Human Plasma	0.1-60 ng/mL	K ₃ EDTA, Lithium Heparin
Aprepitant	LC-MS/MS	Human Plasma	20-6000 ng/mL	EDTA
Atenolol	LC-MS/MS	Human Plasma	5-1500 ng/mL	EDTA
Atomoxetine	LC-MS/MS	Human Plasma	5.2-1560 ng/mL	Heparin
Atorvastatin; Parahydroxy Atorvastatin; Orthohydroxy Atorvastatin	LC-MS/MS	Human Plasma	0.2-80 ng/mL; 0.1-40 ng/mL; 0.02-8 ng/mL	K ₂ and K ₃ EDTA
Atropine	LC-MS/MS	Human Plasma	0.1-20 ng/mL	Heparin
Baclofen	LC-MS/MS	Human Plasma	2-500 ng/mL	EDTA
Baclofen	LC-MS/MS	Mouse Plasma	20-5000 ng/mL	Heparin
Bicalutamide	LC-MS/MS	Human Plasma	10-5000 ng/mL	K ₂ and K ₃ EDTA
Buprenorphine	LC-MS/MS	Human Plasma	5-1000 pg/mL	Heparin
Buprenorphine	LC-MS/MS	Human Plasma	25-5000 pg/mL	Heparin
Buprenorphine; Norbuprenorphine	LC-MS/MS	Human Plasma	10-2500 pg/mL; 10-2500 pg/mL	Sodium Heparin
Buprenorphine; Norbuprenorphine	LC-MS/MS	Human Plasma	5-1000 pg/mL; 10-2000 pg/mL	Heparin
Bupropion	LC-MS/MS	Human Plasma	1-200 ng/mL	EDTA
Bupropion; Hydroxybupropion	LC-MS/MS	Human Plasma	1-400 ng/mL; 4-1600 ng/mL	EDTA
C-Peptide	RIA	Human Serum	0.05-4 ng/mL	N/A
C-Peptide	RIA	Human Serum	100-4000 pg/mL	N/A
Caffeine	LC-MS/MS	Human Plasma	20-5000 ng/mL	Lithium Heparin
Caffeine; Paraxanthine	LC-MS/MS	Human Plasma	20-5000 ng/mL; 20-5000 ng/mL	Lithium Heparin
Candesartan	LC-MS/MS	Human Plasma	0.5-200 ng/mL	EDTA
Candesartan	LC-MS/MS	Human Plasma	2-200 ng/mL	EDTA
Capsaicin	LC-MS/MS	Human Plasma	5-500 pg/mL	Heparin
Celecoxib	LC-MS/MS	Human Plasma	10-2500 ng/mL	EDTA
CEMA; 3-HPMA; HBMA	LC-MS/MS	Human Urine	20-5000 ng/mL; 20-5000 ng/mL; 0.2-150 ng/mL	N/A
Cetirizine	LC-MS/MS	Human Plasma	10-750 ng/mL	EDTA
Cetirizine	LC-MS/MS	Human Plasma	2-400 ng/mL	EDTA
Chlorpheniramine	LC-MS/MS	Human Plasma	0.1-30 ng/mL	Heparin
Cholecalciferol	LC-MS/MS	Human Plasma	0.1-15 ng/mL	EDTA
Clarithromycin; 14-Hydroxy-Clarithromycin	LC-MS/MS	Human Plasma	1-4000 ng/mL; 0.5-2000 ng/mL	EDTA
Clavulanic Acid	LC-MS/MS	Human Plasma	50-5000 ng/mL	EDTA
Clodronate	LC-MS/MS	Human Plasma	3-600 ng/mL	Heparin
Clonazepam	LC-MS/MS	Human Plasma	0.1-20 ng/mL	EDTA
Clopidogrel; Clopidogrel Acid	LC-MS/MS	Human Plasma	0.01-5 ng/mL; 10-5000 ng/mL	EDTA
Codeine	LC-MS/MS	Human Plasma	2-200 ng/mL; 0.1-30 ng/mL	Heparin
Colchicine	LC-MS/MS	Human Plasma	0.03-10 ng/mL	EDTA
Cortisol; 6-Beta-Hydroxycortisol	LC-MS/MS	Human Urine	2-500 ng/mL; 30-2000 ng/mL	N/A

Validated Assay List

Compound Name	Technique	Species/Matrix	Range	Anti-coag.
Cortisol; 6-Beta-Hydroxycortisol	LC-MS/MS	Human Urine	2-250 ng/mL; 25-1000 ng/mL	N/A
Cotinine	LC-MS/MS	Human Saliva	1-100 ng/mL	N/A
Cyclosporine A	LC-MS/MS	Human Whole Blood	5-1000 ng/mL	K ₂ EDTA
Cyproterone Acetate; 15-Beta-Hydroxycyproterone Acetate	LC-MS/MS	Human Plasma	0.1-40 ng/mL	Heparin
Danazol	LC-MS/MS	Human Plasma	0.5-200 ng/mL	K ₃ EDTA
Danazol	LC-MS/MS	Human Serum	0.5-250 ng/mL	N/A
Darifenacin	LC-MS/MS	Human Plasma	0.05-10 ng/mL	EDTA
Deferiprone	LC-MS/MS	Human Serum	0.1-20 ug/mL	N/A
Desipramine	LC-MS/MS	Human Plasma	0.5-200 ng/mL	Heparin
Desmopressin	LC-MS/MS	Human Plasma	2-250 pg/mL	K ₃ EDTA
Dexamethasone	LC-MS/MS	Human Plasma	0.1-50 ng/mL	EDTA
Dexlansoprazole	LC-MS/MS	Human Plasma	20-4000 ng/mL	EDTA
Dextromethorphan; Dextrorphan	LC-MS/MS	Rat Plasma	2-1000 ng/mL; 25-5000 ng/mL	EDTA, Heparin
Dextromethorphan; Dextrorphan	LC-MS/MS	Dog Plasma	2-1000 ng/mL; 25-5000 ng/mL	EDTA
Dextromethorphan; Dextrorphan	LC-MS/MS	Human Plasma	0.2-200 ng/mL; 2.5-2500 ng/mL	Heparin
Dextromethorphan; Dextrorphan	LC-MS/MS	Mouse Plasma	2-1000 ng/mL; 25-5000 ng/mL	EDTA
Dextromethorphan; Dextrorphan (total)	LC-MS/MS	Human Urine	0.12-30 ng/mL; 0.6-300 ng/mL	N/A
Dextromethorphan; Dextrorphan (total); 3-Methoxymorphinan	LC-MS/MS	Human Plasma	0.12-30 ng/mL; 0.6-300 ng/mL; 0.24-30 ng/mL	Lithium Heparin
Diazepam; Nordiazepam	LC-MS/MS	Human Plasma	2-500 ng/mL	EDTA
Diclofenac	LC-MS/MS	Human Plasma	25-5000 pg/mL	K ₂ EDTA
Diclofenac	LC-MS/MS	Human Plasma	0.5-50 ng/mL	Heparin
Diclofenac	LC-MS/MS	Human Plasma	2.5-500 ng/mL	K ₂ EDTA
Digoxin	LC-MS/MS	Human Plasma	10-3500 pg/mL	K ₂ EDTA
Dihydrotestosterone (DHT); Testosterone	LC-MS/MS	Human Serum	10-2000 pg/mL; 50-10000 pg/mL	N/A
Diltiazem; Deacetyldiltiazem; N-Desmethyldiltiazem	LC-MS/MS	Human Plasma	2-100 ng/mL; 0.5-25 ng/mL	Heparin
Diphenhydramine	LC-MS/MS	Human Plasma	0.5-150 ng/mL	EDTA
Dipyridamole	LC-MS/MS	Human Plasma	10-3000 ng/mL	EDTA
Dolasetron; Hydrodolasetron	LC-MS/MS	Human Plasma	1-200 ng/mL	EDTA
Donepezil	LC-MS/MS	Human Plasma	0.3-45 ng/mL	Heparin
Doxylamine	LC-MS/MS	Human Plasma	0.2-100 ng/mL	K ₂ EDTA
Dronedarone	LC-MS/MS	Human Plasma	1-250 ng/mL	EDTA
Drospirenone	LC-MS/MS	Human Plasma	0.5-50 ng/mL	EDTA
Duloxetine	LC-MS/MS	Human Plasma	0.1-20 ng/mL	Heparin
Efavirenz	LC-MS/MS	Human Plasma	10-5000 ng/mL	EDTA
Emtricitabine	LC-MS/MS	Human Plasma	10-2500 ng/mL	Heparin
Enalapril; Enalaprilat	LC-MS/MS	Human Plasma	0.5-100 ng/mL; 0.25-50 ng/mL	EDTA
Eprosartan	LC-MS/MS	Human Plasma	20-4000 ng/mL	EDTA
Eptifibatide	LC-MS/MS	Human Plasma	20-4000 ng/mL	EDTA
Escitalopram	LC-MS/MS	Human Plasma	0.2-40 ng/mL	K ₂ EDTA
Esomeprazole	LC-MS/MS	Human Plasma	2-2000 ng/mL	EDTA

Validated Assay List

Compound Name	Technique	Species/Matrix	Range	Anti-coag.
Estradiol	LC-MS/MS	Human Serum	2.5-500 pg/mL	N/A
Ethinyl Estradiol	LC-MS/MS	Human Plasma	2-500 pg/mL	EDTA
Exemestane	LC-MS/MS	Human Plasma	0.1-50 ng/mL	K ₃ EDTA
Exenatide	EIA	Human Plasma	20-400 pg/mL	Heparin
Ezetimibe-Free	LC-MS/MS	Human Plasma	0.05-10 ng/mL	EDTA
Ezetimibe-Total	LC-MS/MS	Human Plasma	0.5-100 ng/mL	EDTA
Famotidine	LC-MS/MS	Human Plasma	1-200 ng/mL	EDTA
Febuxostat	LC-MS/MS	Human Plasma	15-3000 ng/mL	Sodium Heparin
Felodipine	LC-MS/MS	Human Plasma	10-2000 pg/mL	EDTA
Fenofibrate; Fenofibric Acid; Reduced Fenofibric Acid	LC-MS/MS	Human Plasma	0.05-10 ng/mL; 10-2000 ng/mL; 0.5-100 ng/mL	EDTA
Fentanyl	LC-MS/MS	Human Serum	10-5000 pg/mL	N/A
Fentanyl	LC-MS/MS	Rabbit Serum	100-5000 pg/mL	N/A
Fentanyl	LC-MS/MS	Rat Plasma	20-7500 pg/mL	EDTA
Fentanyl	LC-MS/MS	Human Plasma	10-2000 pg/mL	EDTA
Fexofenadine	LC-MS/MS	Human Plasma	2-850 ng/mL	Heparin
Fexofenadine	LC-MS/MS	Human Plasma	2-400 ng/mL	Heparin
Fluoxetine; Duloxetine; Paroxetine; Sertraline	LC-MS/MS	Human Plasma	1-500 ng/mL; 0.1-50 ng/mL; 0.2-100 ng/mL; 0.5-250 ng/mL	K ₂ EDTA
Flurbiprofen	LC-MS/MS	Human Plasma	50-30000 ng/mL	EDTA
Fluticasone	LC-MS/MS	Human Plasma	1-150 pg/mL	K ₂ EDTA
Fluvastatin	LC-MS/MS	Human Plasma	1-400 ng/mL	K ₂ EDTA
Fluvoxamine	LC-MS/MS	Human Plasma	0.200-100 ng/mL	EDTA
Furosemide	LC-MS/MS	Human Plasma	0.1-20 ug/mL	EDTA
Furosemide	LC-MS/MS	Human Urine	0.1-20 ug/mL	N/A
Galanthamine	LC-MS/MS	Human Plasma	0.3-120 ng/mL	EDTA
Gemcitabine; dFdU	LC-MS/MS	Human Plasma	100-20000 ng/mL; 247-49300 ng/mL	Heparin
Gestodene	LC-MS/MS	Human Plasma	50-10000 pg/mL	EDTA
Glimepiride	LC-MS/MS	Human Plasma	2-400 ng/mL	EDTA
Glimepiride; Trans hydroxy Glimepiride	LC-MS/MS	Human Plasma	1-200 ng/mL; 0.2-40 ng/mL	EDTA
Glucagon	LC-MS/MS	Human Plasma	100-10000 pg/mL	EDTA
Glucagon	LC-MS/MS	Rat Plasma	1-150 ng/mL	EDTA
Glucose	Enzyme Assay	Human Plasma	60-1800 ug/mL	K ₃ EDTA
Glyburide	LC-MS/MS	Human Plasma	2-400 ng/mL	EDTA
Glyburide; cis-3-Hydroxy Glyburide; trans-4-Hydroxy Glyburide	LC-MS/MS	Human Plasma	0.5-200 ng/mL; 0.1-50 ng/mL	Heparin
Glyburide; cis-3-Hydroxy Glyburide; trans-4-Hydroxy Glyburide	LC-MS/MS	Human Urine	5-2500 ng/mL; 20-10000 ng/mL	N/A
HEMA	LC-MS/MS	Human Urine	0.1-20 ng/mL	N/A
Hydrochlorothiazide	LC-MS/MS	Human Plasma	0.5-150 ng/mL	EDTA
Hydrocodone	LC-MS/MS	Human Plasma	0.1-40 ng/mL	Heparin
Ibandronate	LC-MS/MS	Human Serum	0.2-100 ng/mL	N/A
Ibuprofen	LC-MS/MS	Human Plasma	0.25-50 ug/mL	EDTA

Validated Assay List

Compound Name	Technique	Species/Matrix	Range	Anti-coag.
Ibuprofen	LC-MS/MS	Human Plasma	0.5-50 ug/mL	Lithium Heparin
Indapamide	LC-MS/MS	Human Whole Blood	1-200 ng/mL	EDTA
Irbesartan	LC-MS/MS	Human Plasma	30-6000 ng/mL	EDTA
Isradipine	LC-MS/MS	Human Plasma	50-10000 pg/mL	EDTA
Ketoconazole	LC-MS/MS	Human Plasma	20-4000 ng/mL	EDTA
Ketoconazole	LC-MS/MS	Human Plasma	20-15000 ng/mL	Heparin
Ketoprofen	LC-MS/MS	Human Plasma	0.25-50 ng/mL	EDTA, Heparin
Lansoprazole	LC-MS/MS	Human Plasma	2-1200 ng/mL	Heparin
Lenalidomide	LC-MS/MS	Human Plasma	5-1000 ng/mL	K ₂ EDTA
Letrozole	LC-MS/MS	Human Plasma	0.5-100 ng/mL	EDTA
Levonorgestrel	LC-MS/MS	Human Plasma	50-10000 pg/mL	EDTA
Losartan; Losartan Acid Metabolite	LC-MS/MS	Human Plasma	1-300 ng/mL	Sodium EDTA
Lovastatin; Lovastatin Hydroxy-Acid	LC-MS/MS	Human Plasma	0.1-20 ng/mL; 0.1-10 ng/mL	EDTA
Lovastatin; Lovastatin Hydroxy-Acid	LC-MS/MS	Human Plasma	0.1-20 ng/mL	Heparin
Malondialdehyde	EIA	Human Plasma	20-500 ng/mL	K ₂ EDTA
Mecamylamine	LC-MS/MS	Human Plasma	1-100 ng/mL	EDTA
Mesalamine	LC-MS/MS	Human Plasma	3-1000 ng/mL	K ₂ EDTA
Metaxalone	LC-MS/MS	Human Plasma	25-5000 ng/mL	EDTA
Metformin	LC-MS/MS	Human Plasma	30-6000 ng/mL	EDTA
Methotrexate	LC-MS/MS	Human Plasma	2-400 ng/mL	K ₂ EDTA
Methylphenidate	LC-MS/MS	Dog Plasma	0.25-10 ng/mL	Heparin
Methylphenidate	LC-MS/MS	Human Plasma	0.1-20 ng/mL	EDTA
Methylphenidate-d; Methylphenidate-L	LC-MS/MS	Human Plasma	0.5-50 ng/mL; 0.05-5 ng/mL	EDTA
Methylphenidate-d; Ritalinic Acid-d; Methylphenidate-L; Ritalinic Acid-L	LC-MS/MS	Rat Plasma	0.5-50 ng/mL; 5-500 ng/mL	Potassium Oxalate/ Sodium Fluoride
Methylphenidate-d; Ritalinic Acid-d; Methylphenidate-L; Ritalinic Acid-L	LC-MS/MS	Mouse Plasma	1-50 ng/mL; 5-500 ng/mL	Potassium Oxalate/ Sodium Fluoride
Methylprednisolone	LC-MS/MS	Dog Plasma	0.5-50 ng/mL	EDTA
Metoclopramide	LC-MS/MS	Human Plasma	1-150 ng/mL	Heparin
Midazolam; Hydroxymidazolam	LC-MS/MS	Human Plasma	0.1-40 ng/mL; 0.1-20 ng/mL	Heparin
Misoprostol Acid	LC-MS/MS	Human Plasma	5-500 pg/mL	K ₂ EDTA
Mometasone	LC-MS/MS	Human Plasma	0.25-25 pg/mL	EDTA
Monohydroxy-3-Butenyl Mercapturic Acid (MHBMA)	LC-MS/MS	Human Urine	0.1-20 ng/mL	N/A
Montelukast	LC-MS/MS	Human Plasma	3-800 ng/mL	Heparin
Morphine; Morphine-3-Glucuronide; Morphine-6-Glucuronide	LC-MS/MS	Human Serum	0.5-50 ng/mL; 10-1000 ng/mL; 1-100 ng/mL	N/A
Morphine; Morphine-3-Glucuronide; Morphine-6-Glucuronide	LC-MS/MS	Human Plasma	0.25-50 ng/mL; 10-1000 ng/mL	EDTA

Validated Assay List

Compound Name	Technique	Species/Matrix	Range	Anti-coag.
Morphine; Morphine-6-Glucuronide	LC-MS/MS	Human Plasma	0.25-50 ng/mL; 1-200 ng/mL	K ₃ EDTA
Moxifloxacin	LC-MS/MS	Human Plasma	0.5-50 ng/mL	Heparin
NAB-Total; NAT-Total	LC-MS/MS	Human Urine	2-200 pg/mL; 5-1000 pg/mL	N/A
Naftifine	LC-MS/MS	Human Plasma	1-200 ng/mL	K ₃ EDTA
Nalbuphine	LC-MS/MS	Human Plasma	0.5- 50 ng/mL	EDTA
Naloxone	LC-MS/MS	Mouse Plasma	0.5-100 ng/mL	Heparin
Naloxone-Free	LC-MS/MS	Human Plasma	0.2-40 ng/mL	Heparin
Naloxone-Total	LC-MS/MS	Human Plasma	10-2000 pg/mL	EDTA
Naloxone-3-Glucuronide	LC-MS/MS	Human Plasma	75-15000 pg/mL	Heparin
Naltrexone	LC-MS/MS	Dog Plasma	0.1-50 ng/mL	EDTA
Naltrexone	LC-MS/MS	Human Serum	0.1-50 ng/mL	N/A
Naltrexone; 6-Beta-Naltrexol	LC-MS/MS	Human Plasma	4-500 pg/mL; 4-1600 pg/mL	EDTA
Naltrexone; 6-Beta-Naltrexol	LC-MS/MS	Human Serum	10-2000 pg/mL	N/A
Niacin; Nicotinamide; Nicotinuric Acid	LC-MS/MS	Human Plasma	25-5000 ng/mL; 10-2000 ng/mL; 5-1000 ng/mL	EDTA
Nicotine; Cotinine	LC-MS/MS	Human Serum	0.5-50 ng/mL; 1-200 ng/mL	N/A
Nicotine; Cotinine	LC-MS/MS	Human Plasma	0.2-10 ng/mL; 1-100 ng/mL	EDTA
Nicotine; Cotinine	LC-MS/MS	Human Plasma	1-50 ng/mL; 1-100 ng/mL	EDTA
Nicotine; Cotinine; <i>trans</i> -3'-Hydroxycotinine	LC-MS/MS	Human Plasma	0.2-10 ng/mL; 1-100 ng/mL	EDTA, Heparin
Nicotine; Cotinine; <i>trans</i> -3'-Hydroxycotinine; Nicotine Glucuronide; Cotinine Glucuronide; <i>trans</i> -3'-Hydroxycotinine Glucuronide	LC-MS/MS	Human Urine	10-1000 ng/mL; 10-1000 ng/mL; 10-1000 ng/mL; 10-1000 ng/mL; 20-2000 ng/mL; 50-5000 ng/mL	N/A
Nicotine; Cotinine; <i>trans</i> -3'-Hydroxycotinine; Nicotine Glucuronide; Cotinine Glucuronide; <i>trans</i> -3'-Hydroxycotinine Glucuronide	LC-MS/MS	Human Urine	50-5000 ng/mL; 50-5000 ng/mL; 50-5000 ng/mL; 50-5000 ng/mL; 200-20000 ng/mL; 200-20000 ng/mL	N/A
NNAL-Total	LC-MS/MS	Human Urine	5-2000 pg/mL	N/A
NNN; NNAL	LC-MS/MS	Human Urine	0.2-40 pg/mL; 5-2000 pg/mL	N/A
Nordiazepam	LC-MS/MS	Human Plasma	10-500 ng/mL	EDTA
Norethindrone	LC-MS/MS	Human Plasma	0.05-10 ng/mL	K ₂ EDTA
Norgestimate; Deacetylnorgestimate	LC-MS/MS	Human Plasma	20-2500 pg/mL	EDTA
Norgestrel (Reported as Levonorgestrel)	LC-MS/MS	Human Plasma	0.05-10 ng/mL	EDTA, Heparin
Olmesartan	LC-MS/MS	Human Plasma	2.5-500 ng/mL	EDTA
Olmesartan	LC-MS/MS	Human Plasma	2.5-1500 ng/mL	K ₃ EDTA
Omeprazole	LC-MS/MS	Human Plasma	10-2000 ng/mL	Sodium Heparin
Omeprazole	LC-MS/MS	Human Plasma	2-700 ng/mL	EDTA
Ondansetron	LC-MS/MS	Human Plasma	0.5-100 ng/mL	EDTA
Orlistat	LC-MS/MS	Human Plasma	37.5-75000 pg/mL	EDTA
Oxybutynin [R]; Oxybutynin [S]; Desethyloxybutynin [R]; Desethyloxybutynin [S]	LC-MS/MS	Rat Plasma	1-250 ng/mL; 2.5-500 ng/mL	K ₂ EDTA

Validated Assay List

Compound Name	Technique	Species/Matrix	Range	Anti-coag.
Oxybutynin [R]; Oxybutynin [S]; Desethyloxybutynin [R]; Desethyloxybutynin [S]	LC-MS/MS	Human Plasma	0.05-10 ng/mL	EDTA
Oxybutynin; Desethyloxybutynin	LC-MS/MS	Human Plasma	0.05-50 ng/mL	EDTA, Heparin
Oxybutynin; Desethyloxybutynin	LC-MS/MS	Human Plasma	0.05-10 ng/mL	K ₂ EDTA
Oxycodone	LC-MS/MS	Human Plasma	0.2-50 ng/mL	EDTA, Heparin
Oxycodone; Noroxycodone	LC-MS/MS	Human Plasma	0.5-100 ng/mL; 0.02 - 5.00 ng/mL	Sodium EDTA
Oxymorphone	LC-MS/MS	Human Plasma	0.02-5 ng/mL	EDTA
Paliperidone	LC-MS/MS	Human Plasma	0.1-20 ng/mL	K ₂ EDTA
Pantoprazole	LC-MS/MS	Human Plasma	20-6000 ng/mL	Heparin
Paroxetine	LC-MS/MS	Human Plasma	0.1-20 ng/mL	Heparin
Paroxetine	LC-MS/MS	Human Plasma	0.1- 50 ng/mL	EDTA
PEG-IFN alpha 2a (Pegasys)	EIA	Human Serum	250-6000 pg/mL	N/A
Phentermine	LC-MS/MS	Human Plasma	0.25-50 ng/mL	EDTA
Phenytoin	LC-MS/MS	Human Plasma	20-5000 ng/mL	K ₃ EDTA
Pitavastatin	LC-MS/MS	Human Plasma	0.5-150 ng/mL	K ₂ EDTA
Pramipexole	LC-MS/MS	Human Plasma	10-2000 pg/mL	K ₂ EDTA
Pramipexole	LC-MS/MS	Human Plasma	30-7500 pg/mL	K ₂ EDTA
Pravastatin; 3-Alpha-Hydroxy Pravastatin; Pravastatin Lactone	LC-MS/MS	Human Plasma	0.2-100 ng/mL	Heparin
Pravastatin; 3-Alpha Hydroxy Pravastatin; Pravastatin Lactone	LC-MS/MS	Human Urine	0.2-100 ng/mL	N/A
Prednisone; Prednisolone	LC-MS/MS	Human Plasma	0.4-100 ng/mL; 4-1000 ng/mL	EDTA
Progesterone	LC-MS/MS	Human Plasma	0.05-15 ng/mL	EDTA
Progesterone	LC-MS/MS	Human Serum	0.05-15 ng/mL	N/A
Promethazine	LC-MS/MS	Human Plasma	0.1-20 ng/mL	Heparin
Prostaglandin E2 (PGE-2)	EIA	Human Urine	30-1500 pg/mL (STD curve 10-500 pg/mL, samples @ DF=3)	N/A
Pseudoephedrine	LC-MS/MS	Human Plasma	1-200 ng/mL	EDTA
PTH 1-34	ELISA	Human Plasma	6-150 pg/mL	EDTA
PTH 1-84	IRMA	Human Plasma	4 (anchor point) 10-2000 pg/mL	EDTA
Quetiapine	LC-MS/MS	Human Plasma	1-200 ng/mL	EDTA
Quetiapine	LC-MS/MS	Human Plasma	5-1000 ng/mL	EDTA
Quinidine	HPLC/FL	Human Plasma	0.05-10 ug/mL	Heparin
Quinidine	HPLC/FL	Human Plasma	2-250 ng/mL	Heparin
Rabeprazole	LC-MS/MS	Human Plasma	5-400 ng/mL	Heparin
Raloxifene	LC-MS/MS	Human Plasma	5-1000 pg/mL	Heparin
Ramipril; Ramiprilat	LC-MS/MS	Human Plasma	0.1-20 ng/mL	EDTA
Repaglinide	LC-MS/MS	Human Plasma	0.3-75 ng/mL	EDTA
Ribavirin	LC-MS/MS	Human Serum	5-2000 ng/mL	N/A
Ribavirin	LC-MS/MS	Human Red Blood Cells	0.1-10 ug/mL	EDTA
Ribavirin	LC-MS/MS	Human Plasma	10-2000 ng/mL	EDTA

Validated Assay List

Compound Name	Technique	Species/Matrix	Range	Anti-coag.
Rifampin	LC-MS/MS	Human Plasma	0.1-20 ug/mL	K ₂ EDTA
Risedronate	LC-MS/MS	Human Serum	5-2500 pg/mL	N/A
Ritalinic Acid-d; Ritalinic Acid-L	LC-MS/MS	Human Plasma	5-500 ng/mL	EDTA
Ritonavir	LC-MS/MS	Human Plasma	5-1000 ng/mL	Sodium EDTA
Rivastigmine	LC-MS/MS	Human Plasma	0.1-40 ng/mL	K ₃ EDTA
R-Methadone; S-Methadone	LC-MS/MS	Human Plasma	20-4000 ng/mL	K ₂ EDTA
R-Omeprazole; S-Omeprazole	LC-MS/MS	Human Plasma	0.5-100 ng/mL	EDTA
R-Warfarin; S-Warfarin	LC-MS/MS	Human Plasma	12.5-2500 ng/mL	EDTA
R-Warfarin; S-Warfarin	LC-MS/MS	Human Plasma	12.5-2500 ng/mL	Sodium Heparin
Ropinirole	LC-MS/MS	Human Plasma	0.05-7.5 ng/mL	EDTA
Rosiglitazone	LC-MS/MS	Human Plasma	2.5-500 ng/mL	EDTA
Rosiglitazone; Desmethylrosiglitazone	LC-MS/MS	Human Plasma	1-400 ng/mL	EDTA
Rosuvastatin	LC-MS/MS	Human Plasma	0.2-100 ng/mL	EDTA
Rotigotine	LC-MS/MS	Human Plasma	0.01-2 ng/mL	EDTA
Scopolamine	LC-MS/MS	Human Plasma	5-500 pg/mL	EDTA
Sildenafil; Desmethylsildenafil	LC-MS/MS	Human Plasma	1-400 ng/mL	EDTA
Simvastatin; Simvastatin Hydroxy Acid	LC-MS/MS	Human Plasma	0.1-20 ng/mL	EDTA
Sitagliptin	LC-MS/MS	Human Plasma	2.5-500 ng/mL	EDTA
Solifenacin	LC-MS/MS	Human Plasma	0.1-20 ng/mL	EDTA
S-Phenyl Mercapturic Acid	LC-MS/MS	Human Urine	0.02-20 ng/mL	N/A
Sumatriptan	LC-MS/MS	Human Plasma	0.1-15 ng/mL	K ₂ EDTA
Sumatriptan	LC-MS/MS	Human Plasma	0.5-100 ng/mL	K ₂ EDTA
Tacrolimus	LC-MS/MS	Human Whole Blood	0.5-50 ng/mL	K ₂ EDTA
Tadalafil	LC-MS/MS	Human Plasma	1-400 ng/mL	EDTA
Tamsulosin	LC-MS/MS	Human Plasma	0.3-30 ng/mL	Lithium Heparin
Telaprevir	LC-MS/MS	Human Plasma	0.1-20 ug/mL	K ₂ EDTA
Telmisartan	LC-MS/MS	Human Plasma	0.5-100 ng/mL	EDTA
Tenofovir	LC-MS/MS	Human Plasma	5-750 ng/mL	Heparin
Testosterone; Methyltestosterone	LC-MS/MS	Human Plasma	0.05-10 ng/mL	N/A
Tetrahydrocannabinol (THC); THC-OH; THC-COOH	LC-MS/MS	Human Plasma	0.05-10 ng/mL; 0.05-2 ng/mL; 0.75-10 ng/mL	Heparin
Thalidomide	LC-MS/MS	Human Plasma	20-2000 ng/mL	Heparin
Theophylline	LC-MS/MS	Human Plasma	100-20000 ng/mL	EDTA
TNF alpha	EIA	Human Plasma	1-32 pg/mL	Lithium Heparin
Tolterodine; 5-Hydroxymethyl Tolterodine	LC-MS/MS	Human Plasma	15-5000 pg/mL; 15-5000 pg/mL	EDTA
Topiramate	LC-MS/MS	Human Plasma	10-1000 ng/mL	EDTA
Topiramate	LC-MS/MS	Human Plasma	10-5000 ng/mL	EDTA
Torsemide	LC-MS/MS	Human Plasma	20-10000 ng/mL	EDTA
Triamcinolone Acetonide	LC-MS/MS	Human Plasma	50-4000 pg/mL	EDTA
Valsartan	LC-MS/MS	Human Plasma	15-7000 ng/mL	EDTA

Validated Assay List

Compound Name	Technique	Species/Matrix	Range	Anti-coag.
Valsartan	LC-MS/MS	Human Serum	15-7000 ng/mL	N/A
Verapamil; Norverapamil	LC-MS/MS	Human Plasma	2.5-250 ng/mL	Heparin
Verapamil; Norverapamil	LC-MS/MS	Human Plasma	1.3-250 ng/mL	Heparin
Voriconazole	LC-MS/MS	Human Plasma	5-2000 ng/mL	EDTA
Zidovudine; Abacavir	LC-MS/MS	Human Plasma	10-3000 ng/mL; 10-3000 ng/mL	EDTA
Zileuton; A66193	LC-MS/MS	Human Plasma	0.05-10 ug/mL	EDTA
Zoledronate	LC-MS/MS	Human Serum	0.5-100 ng/mL	N/A
Zoledronate	LC-MS/MS	Human Urine	0.5-100 ng/mL	N/A
Zolpidem	LC-MS/MS	Human Plasma	1-250 ng/mL	EDTA
Zonisamide	LC-MS/MS	Human Serum	30-3000 ng/mL	N/A

Biomarker Assays

Compound Name	Technique	Species/Matrix	Range	Anti-coag.
Active Ghrelin	ELISA	Human Plasma	50-2000 pg/mL	K ₂ EDTA
Active GLP-1	ELISA	Human Plasma	2-100 pM	K ₂ EDTA
C-Peptide	RIA	Human Serum	0.05-4 ng/mL	N/A
IGF-1	ELISA	Human Plasma	0.09-6 ng/mL	K ₂ EDTA
IGFBP-3	ELISA	Human Plasma	0.7-40 ng/mL	K ₂ EDTA
IL-6	ELISA	Human Serum	3.12-300 pg/mL	N/A
PGE-2	ELISA	Human Urine	10-500 pg/mL	N/A
PTH 1-34	ELISA	Human Plasma	6-150 pg/mL	K ₃ EDTA
PTH 1-84	RIA	Human Plasma	10-2000 pg/mL	K ₃ EDTA
Total Adiponectin	ELISA	Human Serum	3.93-250 ng/mL	N/A

Biosimilar Qualified Assays

PK Assays - Measurement of Drug in Biological Matrices	Technique	Species/Matrix	Range	Anti-coag.
Cetuximab	ELISA	Human Serum	0.3-50 ug/mL	N/A
Rituximab	ELISA	Human Serum	1-50 ug/mL	N/A
Peg-Filgrastim (Neulasta)	ELISA	Human Serum	0.2-10 ng/mL	N/A
Exendin-4	ELISA	Human Plasma	20-400 pg/mL	Heparin
Etanercept (Enbrel)	ELISA	Human Serum	20-1000 ng/mL	N/A
Adalimumab (Humira)	ELISA	Human Serum	80-8000 ng/mL	N/A
ADA Assays - Measurement of Anti-drug Antibody Assays	Technique	Species/Matrix	Sensitivity	Anti-coag.
Cetuximab ADA assay	ELISA	Human Serum	10 ng/mL	N/A
Rituximab ADA assay	ELISA	Human Serum	10 ng/mL	N/A
PTH (1-34) ADA assay	ELISA	Human Serum	206 ng/mL	N/A

Celerion Solution

Over 40 years of experience

conducting First-in-Human, clinical Proof-of-Concept, cardiac safety services (TQT, robust QT), ADME and NDA-enabling clinical pharmacology studies

Global presence with 750 beds

at four global Celerion clinical facilities, in addition to a network of audited specialty sites in North America, Europe, South Korea and Singapore. Celerion scans the globe to find solutions to clients' early clinical research needs

Extensive expertise

on modeling and simulation, study design, medical writing (protocols and reports), clinical data sciences, biostatistics, and PK/PD analysis

Highly qualified scientists

assist in designing robust, efficient clinical studies as well as analyzing and interpreting data both during and post clinical conduct

Bioanalytical services

at two global locations for method development and validation for small and large molecule assays for the analysis of drugs, metabolites, and biomarkers in a broad range of biological fluids

Drug Development Services team

provides regulatory, drug development and program management support to complement Celerion's service offering

